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### Performance Characteristics

**METHOD COMPARISON:** Clinical evaluation was conducted comparing the results obtained using the AZOG One-Step HIV 1/2 Test Device (the AZOG test was the precursor to the current SMI test) to FDA approved HIV-1/2 ELISA assay test kits. The study included 94 specimens: both assays identified 54 negative, 25 HIV-1 positive and 12 HIV-2 positive results. The results demonstrated 100% overall agreement (for a percent concordance of  $\geq 99\%$ ) of the AZOG One-Step HIV 1/2 Test Device (Serum/Plasma/Whole Blood) when compared to FDA approved HIV-1/2 ELISA Test.

#### Reference HIV-1/2 ELISA Method

		HIV-1 Positive	HIV-2 Positive	Negative
AZOG Method	Positive	25	12	0
	Negative	0	0	54

#### SENSITIVITY AND SPECIFICITY:

The AZOG One-Step HIV 1/2 Test Device (precursor to the SMI test) demonstrated a sensitivity of 100% on HIV-1 and HIV-2 samples and a specificity of 100%.

### Bibliography

1. Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. *The origin of HIV-1 isolate HTLV-IIIb*. Nature (1993)3/363:466-9
2. Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ, Albert, J, Fenyo, EM, Biberfeld, G, Zagury, JF and Laure, F. *New human and simian HIV-related retroviruses possess functional transactivator (tat) gene*. Nature (1987) 328:548-550
3. Caetano JA *Immunologic aspects of HIV infection*. Acta Med Port (1991) 4 Suppl 1:52S-58S
4. Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. *New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes*. JAMA (1998) 280(1): 42-48
5. Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. *Natural protection against HIV-1 infection provided by HIV-2*. Science (1995) 268:1612-1615
6. Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. *HIV-2 and natural protection against HIV-1 infection*. Science (1996) 272:1959-1960.